

Lacutamab reference antibody



Catalog Number: 176869

Product name

Lacutamab reference antibody

Specificity

Human KIR3DL2 / CD158k

Antibody description

Anti-KIR3DL2 / CD158k Antibody (Lacutamab)

Preparation

Recombinant expression and purified from CHO cells.

Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal


Ig Type

IgG1

Applications


ELISA

Validations



SDS-PAGE

Anti-KIR3DL2 / CD158k Antibody (lacutamab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



SEC-HPLC

The purity of Anti-KIR3DL2 / CD158k Antibody (lacutamab) is 99%, determined by SEC-HPLC.