

EBI3 antibody

Catalog Number: 110169

Product name

EBI3 antibody

Specificity

Human; other species not tested.

Antibody description

EBI3 Rabbit Polyclonal antibody. Positive WB detected in MCF7 cells. Positive IP detected in MCF-7 cells. Observed molecular weight by Western-blot: 25-35 kDa

Preparation

This antibody was obtained by immunization of EBI3 recombinant protein (Accession Number: NM_005755). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:110169(EBI3 antibody) at dilution of 1:300



IP Result of anti-EBI3 (IP:Catalog No:110169, 4ug; Detection:Catalog No:110169 1:300) with MCF-7 cells lysate 3200ug.