

# Human AK2 / Adenylate kinase 2 (His & GST Tag) recombinant protein



Catalog Number: 501169

## General Information

### Gene Name Synonym

ATP-AMP transphosphorylase 2; ATP:AMP phosphotransferase; Adenylate monophosphate kinase; Adenylate kinase 2, mitochondrial, N-terminally processed

### Protein Construction

A DNA sequence encoding the human AK2 isoform 1 (P54819-1) (Met 1-Ile 239) was fused with the N-terminal polyhistidine-tagged GST tag at the N-terminus.

### Organism

Human

### Expression Host

Baculovirus-Insect Cells

## QC Testing

### Activity

Kinase activity untested

### Purity

> 95 % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per  $\mu\text{g}$  of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at  $-70^{\circ}\text{C}$

### Predicted N terminal

Met

### Molecular Mass

The recombinant human AK2/GST chimera consists of 476 amino acids and has a calculated molecular mass of 54.3 kDa. It migrates as an approximately 50 kDa band in SDS-PAGE under reducing conditions.

### Formulation

Supplied as sterile 20mM Tris, 500mM NaCl, pH 8.0, 10% gly, 3mM DTT

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$ . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Human AK2 / Adenylate kinase 2 Protein (His & GST Tag) SDS-PAGE