

# SGLT2 antibody

Catalog Number: 115169

## Product name

SGLT2 antibody

## Specificity

Human; other species not tested.

## Antibody description

SGLT2 Rabbit Polyclonal antibody. Positive WB detected in A549 cells. Positive IHC detected in human kidney tissue. Observed molecular weight by Western-blot: 70 kDa

## Preparation

This antibody was obtained by immunization of SGLT2 recombinant protein (Accession Number: XM\_024450402). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


## Validations



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:115169(SLC5A2 Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human kidney tissue slide using Catalog No:115169(SLC5A2 Antibody) at dilution of 1:200 (under 10x lens)



Immunohistochemistry of paraffin-embedded human kidney tissue slide using Catalog No:115169(SLC5A2 Antibody) at dilution of 1:200 (under 40x lens)