

Anti-FH/Fumarase (2A3) Mouse antibody



Catalog Number: 168169

Product name

WB, ICC/IF

Anti-FH/Fumarase (2A3) Mouse antibody

Dilutions

Specificity

WB: 1/1000-2000

Human, Mouse, Rat

Antibody description

FH/Fumarase (2A3) Mouse monoclonal antibody

Preparation

Antigen: Recombinant Protein of FH Fumarase

Formulation

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality


Monoclonal

Ig Type


Mouse IgG1

Applications

Validations



Western blot analysis of 1) Mouse Brain Tissue, 2) Rat Brain tissue, diluted at 1:2000.



Western blot detection of Fumarase in K562, Ramos, Jurkat, HeLa cell lysates using Fumarase (2A3) Mouse mAb(1:1000 diluted). Predicted band size: 50KDa. Observed band size: 50KDa.