

# CDKAL1 antibody



Catalog Number: 109169

## Product name

CDKAL1 antibody

## Immunogen

[Human CDKAL1 Recombinant protein \(His tag\)](#)

## Specificity

Human, Mouse; other species not tested.

## Antibody description

CDKAL1 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in mouse lung tissue, HepG2 cells, HT-1080 cells, Jurkat cells. Observed molecular weight by Western-blot: 60-65 kDa

## Preparation

This antibody was obtained by immunization of CDKAL1 recombinant protein (Accession Number: NM\_017774). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IF: 1:20-1:200

## Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:109169(CDKAL1 antibody) at dilution of



Immunofluorescent analysis of HepG2 cells using Catalog No:109169(CDKAL1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)