

FCN1 antibody



Catalog Number: 110569

Product name

FCN1 antibody

Specificity

Human; other species not tested.

Antibody description

FCN1 Rabbit Polyclonal antibody. Positive WB detected in human liver tissue, K-562 cells. Observed molecular weight by Western-blot: 35 kDa

Preparation

This antibody was obtained by immunization of FCN1 recombinant protein (Accession Number: NM_002003). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:110569(FCN1 antibody) at dilution of 1:500