

Catalog Number: 101569

## Product name

Anti-Biotinidase / biotinase / BTD antibody

## Immunogen

[Human Biotinidase / biotinase / BTD \(His Tag\) recombinant protein](#)

## Specificity

Human Biotinidase / biotinase / BTD

## Antibody description

Rabbit monoclonal to Biotinidase / biotinase / BTD

## Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Human Biotinidase / biotinase / BTD (rh Biotinidase / biotinase / BTD; P43251; Met1-Asp543).

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC-P


## Dilutions

ELISA: 0.1-0.2 µg/mL

This antibody can be used at 0.1-0.2 µg/mL µg/ml with the appropriate secondary reagents to detect Human Biotinidase / biotinase / BTD.


IHC-P: 1-10 µg/mL

## Validations



Biotinidase / biotinase / BTD Antibody, Rabbit MAb, Immunohistochemistry

Immunochemical staining of human BTD in human liver with rabbit monoclonal antibody (5 µg/mL, formalin-fixed paraffin embedded sections).



Biotinidase / biotinase / BTD Antibody, Rabbit MAb, Immunohistochemistry

Immunochemical staining of human BTD in human kidney with rabbit monoclonal antibody (5 µg/mL,

# Anti-Biotinidase / biotinase / BTD antibody



Catalog Number: 101569

---

formalin-fixed paraffin embedded sections).