

PANK1 antibody



Catalog Number: 113569

Product name

PANK1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PANK1 Rabbit Polyclonal antibody. Positive WB detected in Raji cells, HeLa cells, human heart tissue, MCF7 cells, SW 1990 cells. Observed molecular weight by Western-blot: 64kd

Preparation

This antibody was obtained by immunization of PANK1 recombinant protein (Accession Number: XM_017016338). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:113569(PANK1 antibody) at dilution of 1:300