

R-cadherin antibody



Catalog Number: 114569

Product name

R-cadherin antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

R-cadherin Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, HepG2 cells. Observed molecular weight by Western-blot: 100-120 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001794). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114569(CDH4 antibody) at dilution of 1:300