

Anti-BDKRB2 Rabbit antibody



Catalog Number: 167569

Product name

Anti-BDKRB2 Rabbit antibody

Specificity

Human, Mouse

Antibody description

BDKRB2 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human BDKRB2

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


WB, IP

Dilutions

WB: 1/1000-1/5000

IP: 1/20-1/50

Validations



Western blot detection of BDKRB2 in A549,HL-60,U2OS cell lysates using BDKRB2 Rabbit mAb(1:1000 diluted).Predicted band size:45kDa.Observed band size:80-100kDa.