

Anti-Phospho-PKA RIIalpha(Ser99) Rabbit antibody



Catalog Number: 168569

Product name

Anti-Phospho-PKA RIIalpha(Ser99) Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Phospho-PKA RIIalpha(Ser99) Rabbit monoclonal antibody

Preparation

Antigen: A synthetic phosphopeptide corresponding to residues surrounding Ser99 of human PKA R2

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, IP

Dilutions


1:1000

1:100

1:100

1:20

Validations



Western blot detection of Phospho-PKA RIIalpha (Ser99) in Rat Brain, C6 cell lysates using Phospho-PKA RIIalpha (Ser99) Rabbit mAb(1:1000 diluted). Predicted band size: 46kDa. Observed band size: 50kDa.