

Giantin antibody

Catalog Number: 110969

Product name

Giantin antibody

Specificity

Human; other species not tested.

Antibody description

Giantin Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells, HeLa cells.

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_004487).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, IF

Dilutions

Recommended Dilution:

IF: 1:20-1:200

Validations

Immunofluorescent analysis of HepG2 cells using Catalog No:110969(Giantin Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)