

PLK2 antibody

Catalog Number: 113969

Product name

PLK2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PLK2 Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 75-80 kDa

Preparation

This antibody was obtained by immunization of PLK2 recombinant protein (Accession Number: NM_006622). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113969(PLK2 antibody) at dilution of 1:1000



IP Result of anti-PLK2 (IP:Catalog No:113969, 4ug; Detection:Catalog No:113969 1:300) with HeLa cells lysate 2480ug.