

# Anti-AKR1B1 antibody



Catalog Number: 105969

## Product name

Monoclonal

Anti-AKR1B1 antibody

## Ig Type

## Immunogen

Rabbit IgG

[Human AKR1B1 \(His Tag\) recombinant protein](#)

## Applications

## Specificity

FCM

Human AKR1B1

## Dilutions

## Antibody description

FCM: 0.5-2 µg/Test

Rabbit monoclonal to AKR1B1

## Validations

## Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Human AKR1B1 (rh AKR1B1; P15121; Met1-Phe316).

## Formulation

0.2 µm filtered solution in PBS


## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality



AKR1B1 Antibody, Rabbit MAb, Flow Cytometry

Flow cytometric analysis of Human AKR1B1 expression on HeLa cells. The cells were treated according to manufacturer's manual (BD Pharmingen™ Cat. No. 554714), stained with purified anti-Human AKR1B1, then a FITC-conjugated second step antibody. The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of intact cells.