

# SRA1 antibody



Catalog Number: 106969

---

## Product name

SRA1 antibody

## Specificity

Human; other species not tested.

## Antibody description

SRA1 Rabbit Polyclonal antibody. Positive WB detected in A431 cells, MCF7 cells. Observed molecular weight by Western-blot: 37kd

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: BC156760). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:115661(SRA1 antibody) at dilution of 1:600