

TMOD1 antibody



Catalog Number: 116269

Product name

TMOD1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TMOD1 Rabbit Polyclonal antibody. Positive WB detected in mouse skeletal muscle tissue, human brain tissue, mouse brain tissue, mouse heart tissue. Observed molecular weight by Western-blot: 41kd

Preparation

This antibody was obtained by immunization of TMOD1 recombinant protein (Accession Number: NM_001166116). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:116269(TM0D1 Antibody) at dilution of 1:600