

# ZW10 antibody

Catalog Number: 117269

## Product name

ZW10 antibody

## Specificity

Human; other species not tested.

## Antibody description

ZW10 Rabbit Polyclonal antibody. Positive IP detected in Jurkat cells. Positive WB detected in Jurkat cells, HeLa cells, K-562 cells. Observed molecular weight by Western-blot: 88-90 kDa

## Preparation

This antibody was obtained by immunization of ZW10 recombinant protein (Accession Number: BC128264). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:117269(ZW10 Antibody) at dilution of 1:1000



IP Result of anti-ZW10 (IP:Catalog No:117269, 4ug; Detection:Catalog No:117269 1:500) with Jurkat cells lysate 3200ug.