

ASPRV1 antibody



Catalog Number: 108269

Product name

ASPRV1 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

ASPRV1 Rabbit Polyclonal antibody. Positive WB detected in mouse skin tissue. Observed molecular weight by Western-blot: 28kd,37kd

Preparation

This antibody was obtained by immunization of ASPRV1 recombinant protein (Accession Number: BC031997). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse skin tissue were subjected to SDS PAGE followed by western blot with Catalog No:108269(ASPRV1 antibody) at dilution of 1:300