

PDE7B antibody

Catalog Number: 113669

Product name

PDE7B antibody

IHC: 1:20-1:200

IF: 1:10-1:100


Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PDE7B Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, mouse brain tissue. Positive IHC detected in human brain tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 50kd

Validations



Preparation

This antibody was obtained by immunization of PDE7B recombinant protein (Accession Number: NM_018945). Purification method: Antigen affinity purified.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113669(PDE7B antibody) at dilution of 1:1000

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC, IF


Dilutions

Recommended Dilution:

WB: 1:500-1:5000



Immunohistochemical of paraffin-embedded human brain using Catalog No:113669(PDE7B antibody) at dilution of 1:100 (under 10x lens)




Immunohistochemical of paraffin-embedded human brain using Catalog No:113669(PDE7B antibody) at dilution of 1:100 (under 10x lens)

PDE7B antibody

Catalog Number: 113669

antibody) at dilution of 1:100 (under 40x lens)



Immunofluorescent analysis of HepG2 cells, using PDE7B antibody Catalog No:113669 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).