

RDH13 antibody

Catalog Number: 114669

Product name

RDH13 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RDH13 Rabbit Polyclonal antibody. Positive WB detected in human heart tissue, Y79 cells. Observed molecular weight by Western-blot: 36kd

Preparation

This antibody was obtained by immunization of RDH13 recombinant protein (Accession Number: NM_138412). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:114669(RDH13 antibody) at dilution of 1:500