

# TSPYL4 antibody

Catalog Number: 116369

## Product name

TSPYL4 antibody

## Specificity

Mouse,Rat; other species not tested.

## Antibody description

TSPYL4 Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells.

## Preparation

This antibody was obtained by immunization of TSPYL4 recombinant protein (Accession Number: BC009116). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116369(TSPYL4 Antibody) at dilution of 1:1000