

Human ATG3/APG3/APG3L recombinant protein



Catalog Number: 518769

General Information

Gene Name Synonym

Ubiquitin-like-conjugating enzyme ATG3;
Autophagy-related protein 3; APG3-like; hApg3;
Protein PC3-96; ATG3; APG3; APG3L

Protein Construction

A DNA sequence encoding the Human
ATG3/APG3/APG3L Met1-Met314 is expressed.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μ g protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 μ m filtered solution of 20mM
TrisHCl, 150mM NaCl, 10% Glycerol, pH 8.0.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human ATG3/APG3/APG3L recombinant protein