

# DNA2 antibody

Catalog Number: 110079

## Product name

DNA2 antibody

IF: 1:10-1:100


## Specificity

Human, Mouse; other species not tested.

## Antibody description

DNA2 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive IHC detected in human cervical cancer tissue. Positive WB detected in mouse thymus tissue. Observed molecular weight by Western-blot: 127 kDa

## Validations




## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_001080449). Purification method: Antigen affinity purified.

mouse thymus tissue were subjected to SDS PAGE followed by western blot with Catalog No:110079(DNA2 antibody) at dilution of 1:500

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human cervical cancer using Catalog No:110079(DNA2 antibody) at dilution of 1:50 (under 10x lens)

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC, IF, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

IHC: 1:20-1:200




Immunohistochemical of paraffin-embedded human cervical cancer using Catalog

# DNA2 antibody

Catalog Number: 110079

---

No:110079(DNA2 antibody) at dilution of 1:50  
(under 40x lens)



Immunofluorescent analysis of HepG2 cells using  
Catalog No:110079(DNA2 Antibody) at dilution of  
1:25 and Rhodamine-Goat anti-Rabbit IgG