

# GLUT4 antibody

Catalog Number: 111079

## Product name

GLUT4 antibody

## Specificity

Human; other species not tested.

## Antibody description

GLUT4 Rabbit Polyclonal antibody. Positive IHC detected in human skeletal muscle tissue, human heart tissue.

## Preparation

This antibody was obtained by immunization of GLUT4 recombinant protein (Accession Number: NM\_001042). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC

## Dilutions


Recommended Dilution:

IHC: 1:20-1:200

## Validations



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue slide using Catalog No:111079 (GLUT4 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue slide using Catalog No:111079 (GLUT4 Antibody) at dilution of 1:50 (under 40x lens)