

NDUFC1 antibody

Catalog Number: 113079

Product name

NDUFC1 antibody

Specificity

Human; other species not tested.

Antibody description

NDUFC1 Rabbit Polyclonal antibody. Positive IHC detected in human heart tissue, human kidney tissue.

Preparation

This antibody was obtained by immunization of NDUFC1 recombinant protein (Accession Number: NM_001184986). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human heart tissue slide using Catalog No:113079(NDUFC1 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human heart tissue slide using Catalog No:113079(NDUFC1 Antibody) at dilution of 1:50 (under 40x lens)