

# TM7SF3 antibody



Catalog Number: 116079

## Product name

TM7SF3 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

TM7SF3 Rabbit Polyclonal antibody. Positive WB detected in mouse kidney tissue, mouse pancreas tissue, mouse testis tissue. Observed molecular weight by Western-blot: 64kd

## Preparation

This antibody was obtained by immunization of TM7SF3 recombinant protein (Accession Number: NM\_016551). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

## Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:116079(TM7SF3 antibody) at dilution of 1:1000