

Anti-ENO1 (3F7) Mouse antibody



Catalog Number: 168079

Product name

Anti-ENO1 (3F7) Mouse antibody

Ig Type

Mouse IgG1

Specificity

Human, Mouse, Rat, Monkey

Applications

WB

Antibody description

ENO1 (3F7) Mouse monoclonal antibody

Dilutions

WB: 1/2000

Preparation

Antigen: Purified recombinant human Enolase-1 protein fragments expressed in E.coli.

Formulation

PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.


Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Validations



Western blot detection of Enolase-1 in T47D, 3T3, COS7, C6 and Hela cell lysates using Enolase-1 mouse mAb(dilution 1:2000). Predicted band size: 47kDa. Observed band size: 47kDa.