

# Anti-ENO1 (3F7) Mouse antibody



Catalog Number: 168079

## Product name

Anti-ENO1 (3F7) Mouse antibody

## Specificity

Human, Mouse, Rat, Monkey

## Antibody description

ENO1 (3F7) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human Enolase-1 protein fragments expressed in E.coli.

## Formulation

PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1


## Applications

WB

## Dilutions

WB: 1/2000

## Validations



Western blot detection of Enolase-1 in T47D,3T3,COS7,C6 and HeLa cell lysates using Enolase-1 mouse mAb(dilution 1:2000).Predicted band size:47kDa.Observed band size:47kDa.