

ANKRD53 antibody



Catalog Number: 108079

Product name

ANKRD53 antibody

IHC: 1:20-1:200

IF: 1:50-1:500

Specificity

Human; other species not tested.

Antibody description

ANKRD53 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HepG2 cells. Positive IF detected in HepG2 cells. Positive IHC detected in human kidney tissue, human testis tissue. Observed molecular weight by Western-blot: 65-70 kDa

Preparation

This antibody was obtained by immunization of ANKRD53 recombinant protein (Accession Number: NM_024933). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IF

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:108079(ANKRD53 Antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:108079(ANKRD53 Antibody) at dilution of 1:50 (under 40x lens)




ANKRD53 antibody



Catalog Number: 108079

Immunohistochemistry of paraffin-embedded human kidney tissue slide using Catalog No:108079(ANKRD53 Antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of (-20oc Ethanol) fixed HepG2 cells using Catalog No:108079(ANKRD53 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)