

RASGRP3 antibody



Catalog Number: 114479

Product name

RASGRP3 antibody

Specificity

Human; other species not tested.

Antibody description

RASGRP3 Rabbit Polyclonal antibody. Positive IHC detected in human ovary tumor tissue. Positive WB detected in RAW264.7 cells, PC-3 cells. Observed molecular weight by Western-blot: 70 kDa

Preparation

This antibody was obtained by immunization of RASGRP3 recombinant protein (Accession Number: XM_011532748). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations




RAW264.7 cells were subjected to SDS PAGE followed by western blot with Catalog No:114479(RASGRP3 antibody) at dilution of



1:800

Immunohistochemistry of paraffin-embedded human ovary tumor slide using Catalog No:114479(RASGRP3 Antibody) at dilution of 1:50



Immunohistochemistry of paraffin-embedded human ovary tumor slide using Catalog No:114479(RASGRP3 Antibody) at dilution of 1:50