

Anti-KLRC4-KLRK1 antibody



Catalog Number: 176479

Product name

Anti-KLRC4-KLRK1 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to KLRC4-KLRK1

Preparation

This antigen of this antibody was synthetic peptide from human protein at aa range: 167-216

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

IHC-P

Dilutions

IHC-P: 1:50-200

Validations




Fig1.; Immunohistochemical analysis of paraffin-embedded human-kidney, antibody was diluted at 1:200