

Anti-NTPCR Rabbit antibody



Catalog Number: 167479

Product name

Anti-NTPCR Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

NTPCR Rabbit polyclonal antibody

Preparation

Antigen: Recombinant protein of human Nucleoside triphosphate phosphohydrolase

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, IP


Dilutions

WB: 1/1000

ICC/IF: 1/20

IP: 1/20

Validations



Western blot detection of Nucleoside triphosphate phosphohydrolase in K562,C6,3T3,Hela cell lysates using Nucleoside triphosphate phosphohydrolase Rabbit pAb(1:1000 diluted).Predicted band size:21kDa.Observed band size:21kDa.