

Beta adducin antibody

Catalog Number: 108479

Product name

Beta adducin antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Beta adducin Rabbit Polyclonal antibody. Positive WB detected in K-562 cells. Positive IP detected in K-562 cells. Observed molecular weight by Western-blot: 110 kDa

Preparation

This antibody was obtained by immunization of Beta adducin recombinant protein (Accession Number: NM_001185054). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:117124(ADD2 antibody) at dilution of 1:1000



IP Result of anti-ADD2 (IP:Catalog No:117124, 4ug; Detection:Catalog No:117124 1:1000) with K-562 cells lysate 3200ug.