

Anti-CDH26 antibody

Catalog Number: 175879

Product name

Anti-CDH26 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to CDH26

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human cdh26. aa range:11-60

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IF

Dilutions

WB: 1:500-1:2000

IF: 1:100-1:500

Validations




Fig1.; Western Blot analysis of various cells using Cadherin-26 Polyclonal Antibody