

Anti-LXR alpha Rabbit antibody



Catalog Number: 166879

Product name

Anti-LXR alpha Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

LXR alpha Rabbit polyclonal antibody

Preparation

Antigen: A synthesized peptide derived from human LXR alpha

Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC


Dilutions

WB: 1/500-1/2000

ICC/IF: 1/50-1/200

FC: 1/50

Validations



Western blot analysis of LXR alpha expression in Jurkat cell lysate.