

# XYLT1 antibody

Catalog Number: 116879

---

## Product name

XYLT1 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

XYLT1 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, BxPC-3 cells. Observed molecular weight by Western-blot: 120kd

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_022166). Purification method: Antigen affinity purified.

## Formulation

0.1M NaHCO<sub>3</sub>, 0.1M glycine, 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116879(XYLT1 antibody) at dilution of 1:300