

CSK antibody



Catalog Number: 107179

Product name

CSK antibody

Specificity

Human; other species not tested.

Antibody description

CSK Mouse Monoclonal antibody. Positive IF detected in A431 cells. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 50 kDa

Preparation

This antibody was obtained by immunization of CSK recombinant protein (Accession Number: NM_001354988). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:107179(CSK antibody) at dilution of 1:1000



Immunofluorescent analysis of A431 cells using Catalog No:107179(CSK Antibody) at dilution of 1:50 and Rhodamine-Goat anti-Mouse IgG