

OPRM1 Polyclonal Antibody



Catalog Number: 166579

Product name

OPRM1 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to OPRM1

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 300 to the C-terminus of human OPRM1 (NP_000905.3).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200

Validations



Immunohistochemistry - OPRM1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human placenta using OPRM1 antibody at dilution of 1:100 (40x lens).