

ARID2 Polyclonal Antibody

Catalog Number: 165979

Product name

ARID2 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to ARID2

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1586-1835 of human ARID2 (NP_689854.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions


WB 1:500 - 1:2000

Validations



Western blot - ARID2 Polyclonal Antibody

Western blot analysis of extracts of LO2 cells, using ARID2 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL



Basic Kit .Exposure time: 30s.

Western blot - ARID2 Polyclonal Antibody

Western blot analysis of extracts from normal (control) and ARID2 knockout HeLa cells, using ARID2 antibody at 1:500 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit .Exposure time: 10s.