

# Anti-Phospho-RSK4 (Ser232) Rabbit antibody



Catalog Number: 168979

## Product name

Anti-Phospho-RSK4 (Ser232) Rabbit antibody

## Specificity

Human, Mouse

## Antibody description

Phospho-RSK4 (Ser232) Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic phosphopeptide corresponding to residues surrounding Ser235 of human RSK4

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, ICC/IF


## Dilutions

1:2000

1:20

1:50

## Validations



Western blot detection of Phospho-RSK4 (Ser232) in HeLa cell lysates using Phospho-RSK4 (Ser232) Rabbit mAb(1:1000 diluted). Predicted band size:84kDa.Observed band size:84kDa.