

H12N1 HA (His Tag) recombinant protein



Catalog Number: 502279

General Information

Protein Construction

A DNA sequence encoding the extracellular domain of influenza A virus hemagglutinin (A/mallard duck/Alberta/342/1983 (H12N1)) (ABB88099.1) (Met 1-Lys 527) was fused with a polyhistidine tag at the C-terminus.

Organism

H12N1

Expression Host

Human Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Asp 18

Molecular Mass

The secreted recombinant hemagglutinin of influenza A H12N1 HA (A/mallard duck/Alberta/342/1983(H12N1)) comprises 521

amino acids and has a predicted molecular mass of 59.3 kDa. As a result of glycosylation, it migrates as an approximately 85 kDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Influenza A H12N1 (A/mallard duck/Alberta/342/1983) Hemagglutinin / HA Protein (His Tag) SDS-PAGE