

# PRRX2 antibody

Catalog Number: 114279

## Product name

PRRX2 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

PRRX2 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 27-34 kDa

## Preparation

This antibody was obtained by immunization of PRRX2 recombinant protein (Accession Number: NM\_016307). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114279(PRRX2 Antibody) at dilution of 1:600



IP Result of anti-PRRX2 (IP:Catalog No:114279, 4ug; Detection:Catalog No:114279 1:500) with mouse brain tissue lysate 3440ug.