

# Anti-Monoamine Oxidase A/MAOA Rabbit antibody



Catalog Number: 167279

## Product name

Rabbit IgG

Anti-Monoamine Oxidase A/MAOA Rabbit antibody

## Applications

## Specificity

WB, IHC-P, ICC/IF, FC

Human, Mouse, Rat

## Dilutions

## Antibody description

WB: 1/1000

Monoamine Oxidase A/MAOA Rabbit polyclonal antibody

IHC: 1/20

ICC/IF: 1/20

## Preparation

FC: 1/20

Antigen: A synthetic peptide of human Monoamine Oxidase A

## Validations

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.


## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type



Western blot detection of Monoamine Oxidase A in C6 cell lysates using Monoamine Oxidase A Rabbit pAb(1:1000 diluted). Predicted band size:60kDa. Observed band size:60kDa.