

FKBP6 antibody



Catalog Number: 110679

Product name

FKBP6 antibody

Specificity

Human; other species not tested.

Antibody description

FKBP6 Rabbit Polyclonal antibody. Positive IF detected in MCF-7 cells. Positive WB detected in human testis tissue, HEK-293 cells. Observed molecular weight by Western-blot: 36kd

Preparation

This antibody was obtained by immunization of FKBP6 recombinant protein (Accession Number: NM_003602). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:110679(FKBP6 antibody) at dilution of 1:600



Immunofluorescent analysis of MCF-7 cells, using FKBP6 antibody Catalog No:110679 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent DNA dye).