

Anti-CD6 / Cluster of Differentiation 6 antibody



Catalog Number: 101679

Product name

Anti-CD6 / Cluster of Differentiation 6 antibody

Immunogen

[Mouse CD6 / Cluster of Differentiation 6 \(His Tag\) recombinant protein](#)

Specificity

Mouse CD6 / Cluster of Differentiation 6

Antibody description

Rabbit monoclonal to CD6 / Cluster of Differentiation 6

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Mouse CD6 / Cluster of Differentiation 6 (rM CD6 / Cluster of Differentiation 6; Q91WN5; Met1-Gly396).

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


FCM

Dilutions

IHC-P: 1-10 µg/mL

FCM: 0.5-2 µg/Test

Validations



CD6 / Cluster of Differentiation 6 Antibody, Rabbit MAb, Flow Cytometry

Flow cytometric analysis of CD6 expression on BABL/c splenocytes. Cells were stained with purified anti-Mouse CD6, then a FITC-conjugated second step antibody. The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of intact cells.