

PDHB antibody



Catalog Number: 113679

Product name

PDHB antibody

WB: 1:500-1:5000

IP: 1:200-1:1000

Specificity


Human, Mouse, Rat; other species not tested.

IHC: 1:20-1:200

Antibody description

PDHB Rabbit Polyclonal antibody. Positive IHC detected in human skin cancer tissue. Positive IP detected in mouse heart tissue. Positive WB detected in human heart tissue, human skeletal muscle tissue, K-562 cells, mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 34 kDa

Validations




Preparation

This antibody was obtained by immunization of PDHB recombinant protein (Accession Number: NM_000925). Purification method: Antigen affinity purified.

human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:113679(PDHB antibody) at dilution of 1:300

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human skin cancer using Catalog No:113679(PDHB antibody) at dilution of 1:100 (under 10x lens)

Ig Type


Rabbit IgG

Applications

ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:



PDHB antibody

Catalog Number: 113679

Immunohistochemical of paraffin-embedded human skin cancer using Catalog No:113679(PDHB antibody) at dilution of 1:100 (under 40x lens)



IP Result of anti-PDHB (IP:Catalog No:113679, 4ug; Detection:Catalog No:113679 1:300) with mouse heart tissue lysate 4000ug.