

REEP1 antibody

Catalog Number: 114679

Product name

REEP1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

REEP1 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, rat brain tissue. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 22kd

Preparation

This antibody was obtained by immunization of REEP1 recombinant protein (Accession Number: NM_022912). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114679(REEP1 antibody) at dilution of 1:1000



IP Result of anti-REEP1 (IP:Catalog No:114679, 3ug; Detection:Catalog No:114679 1:500) with mouse brain tissue lysate 4000ug.