

Anti-Lactate Dehydrogenase C Rabbit antibody



Catalog Number: 168679

Product name

Anti-Lactate Dehydrogenase C Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Human

Applications

WB

Antibody description

Lactate Dehydrogenase C Rabbit monoclonal antibody

Dilutions

1:1000-1:5000


Preparation

Antigen: A synthetic peptide of human Lactate Dehydrogenase C

Validations

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.



Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Western blot detection of Lactate Dehydrogenase C in HeLa cell lysates using Lactate Dehydrogenase C Rabbit mAb(1:1000 diluted). Predicted band size:36kDa.Observed band size:36kDa.

Clonality

Monoclonal