

# E-selectin /CD62E antibody



Catalog Number: 110379

---

## Product name

E-selectin /CD62E antibody

## Specificity

Human,Mouse,Rat; other species not tested.

## Antibody description

E-selectin /CD62E Rabbit Polyclonal antibody.  
Positive WB detected in Jurkat cells, L02 cells.  
Observed molecular weight by Western-blot: 67 kDa

## Preparation

This antibody was obtained by immunization of E-selectin /CD62E recombinant protein (Accession Number: NM\_000450). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:110379(SELE antibody) at dilution of 1:300